

AD-A118 424

RESEARCH INST OF NATIONAL DEFENCE STOCKHOLM (SWEDEN) F/G 19/1  
ANALYSES OF EXPLOSIVES. PHYSICAL-CHEMICAL DATA, CHROMATOGRAMS, --ETC(U)  
OCT 78 A ALM, O DALMAN, I FROELEN-LINDGREN  
FOA-C-20267-D1

NL

UNCLASSIFIED

TOP  
SECRET



END

DATE

FILED

68-82

DTIC

NATIONAL DEFENCE RESEARCH INSTITUTE  
Department 2  
S-104 50 Stockholm

FOA-Report  
C-20267-D1  
October 1978

1

# ANALYSES OF EXPLOSIVES

PHYSICAL-CHEMICAL DATA, CHROMATOGRAMS, MASS-, IR- AND NMR-  
SPECTRA. COLOUR REACTIONS AND THIN LAYER CHROMATOGRAPHY

Anita Alm, Olof Dalman, Inger Frölen-Lindgren,  
Felix Hultén, Tom Karlsson and Monica Kowalska


81 pp

## Summary

The present report comprises analytical data of frequently used explosives and tri-, di- and mononitrotoluenes. The data given are obtained by instrumental methods such as MS (70 eV, 20 eV, CI), IR and NMR ( $^1\text{H}$ ,  $^{13}\text{C}$ ), in the form of spectra, and HPLC and GC, as chromatograms. Some data of noninstrumental methods such as colour reactions, thin layer chromatographic applications and pertinent physical-chemical data of these substances are also given.

Uppdragsnummer: D1 14

Accession No.	
100-10001	<input checked="" type="checkbox"/>
100-10002	<input type="checkbox"/>
100-10003	<input type="checkbox"/>
100-10004	<input type="checkbox"/>
FL 288. <i>inf. file</i>	
Distribution/	
Available by Codes	
A and/or	
List Special	
A	



DTIC  
ELECTE  
AUG 23 1982  
A

This document has been approved  
for public release and sale; its  
distribution is unlimited.

<u>CONTENTS</u>	Page
INTRODUCTION	4
PHYSICAL DATA	4
Table 1. List of explosives. Formulas and physical data.	5
Table 2. List of nitrotoluenes. Formulas and physical data.	8
Table 3. Solubilities of explosives.	10
CHROMATOGRAPHY	11
High pressure liquid chromatography (HPLC)	11
Table 4. Retentionvolumes of explosives on $\mu$ -Porasil column	11
Gas chromatography (GC)	16
SPECTROSCOPY	20
Mass spectrometry (MS)	20
Infrared spectroscopy (IR)	20
Table 5. IR-frequencies of (x)-NO <sub>2</sub> associated with nitric esters, nitramines, and nitro compounds.	20
Nuclear magnetic resonance (NMR)	21
SPECTRA	22-71
1. EGDN      Ethylene glycol dinitrate	22
2. NG        Nitroglycerine	24
3. PETN      Pentaerythritol tetranitrate	26
4. NC        Nitrocellulose	28
5. RDX       Hexogen	30
6. HMX       Octogen	32
7. TETR      Tetryl	34
8. AM-PIKR   Ammonium picrate	36
9. HNS       Hexanitrostilbene	38
10. TNB      Trinitrobenzene	40
11. TNT      Trotyl, 2,4,6-Trinitrotoluene	42

12. 2,4,5-TNT	2,4,5-Trinitrotoluene	44
13. 2,3,6-TNT	2,3,6- "	46
14. 2,3,5-TNT	2,3,5- "	48
15. 2,3,4-TNT	2,3,4- "	50
16. 3,4,5-TNT	3,4,5- "	52
17. 2,6-DNT	2,6-Dinitrotoluene	54
18. 2,5-DNT	2,5- "	56
19. 2,4-DNT	2,4- "	58
20. 2,3-DNT	2,3- "	60
21. 3,5-DNT	3,5- "	62
22. 3,4-DNT	3,4- "	64
23. 2-MNT	2-Mononitrotoluene	66
24. 3-MNT	3- "	68
25. 4-MNT	4- "	70
COLOUR REACTIONS		72
Spot tests		72
Table 6. Colours with spot test reagents		74
Detection on silica gel layers		75
Table 7. Colours on silica gel layers		76
THIN LAYER CHROMATOGRAPHY		77
Figure 1. Two dimensional chromatogram of explosives		78
Table 8. $R_f$ -values on silica gel layers		79
REFERENCES		81

## INTRODUCTION

Literature searches<sup>1,2</sup> on identification and analyses of explosives have shown that there is no comprehensive report on analytical methods for this group of substances. (A recently published review is Yinon, J: Analysis of explosives. CRC Crit Rev Anal Chem 7 (1977) 1-35.)

Due to this fact the present work was carried out. The report comprises spectra and chromatograms of 25 substances of which eleven are frequently used explosives and the remaining fourteen are nitrotoluene isomers.

Instrumental methods have been applied. The spectral analyses include mass spectrometry (electron impact, chemical ionization), infrared spectrometry and nuclear magnetic resonance, (<sup>1</sup>H, <sup>13</sup>C). The chromatographic analyses include gas chromatography and high pressure liquid chromatography. Data on instrumentation and experimental conditions are given in each chapter.



As a complement, some pertinent physical-chemical data are presented in tabulated form. Data on qualitative analyses are given, based on noninstrumental methods such as colour reactions and thin layer chromatography.

The collection of data presented will hopefully be of assistance in the qualitative and quantitative analysis of the components in explosive compositions.





## PHYSICAL DATA

Physical data given in tables 1-3 have been collected from different sources: compilations such as Landholt-Börnstein, Beilstein, Handbook of chemistry and physics, Fedoroff<sup>4</sup>, books such as Meyer<sup>5</sup>, Bofors<sup>6</sup> and a few published papers.

Table 1. List of explosives. Formulas and physical data.

No	Explosive	Abbrev	Synonym	Formula	Mol. weight	Melt. point	Deflag. temp.	Remark
1	Ethylene glycol dinitrate	EGDN	(Di)nitroglykol, Glykoldinitrat	$\begin{array}{c} \text{H} & & \text{H} \\   & &   \\ \text{H}-\text{C}-\text{O}-\text{NO}_2 \\   & &   \\ \text{H}-\text{C}-\text{O}-\text{NO}_2 \\   & &   \\ \text{H} & & \text{H} \end{array}$	152,1	-22	217	
2	Nitroglycerine	NG	Glycerin- el glyceroltrinitrat, Trinitro- glycerin	$\begin{array}{c} \text{H} & & \text{H} & & \text{H} \\   & &   & &   \\ \text{H}-\text{C}-\text{O}-\text{NO}_2 \\   & &   & &   \\ \text{H}-\text{C}-\text{O}-\text{NO}_2 \\   & &   & &   \\ \text{H}-\text{C}-\text{O}-\text{NO}_2 \\   & &   & &   \\ \text{H} & & \text{H} & & \text{H} \end{array}$	227,1	13,5 2,8	223-25	Modifications: 13,5 °C triclinic, stable 2,8 °C dipyramidal, rhombic, labile
3	Pentaerythritol tetra- nitrate	PETN	Pentaerytrit- el Pentaery- tritol-tetranitrat, Pentrit, Pentryl, Penta, Nitro- penta (-erytrit)	$\begin{array}{c} \text{NO}_2 \\   \\ \text{O} \\   \\ \text{CH}_2 \\   \\ \text{O}_2\text{N}-\text{O}-\text{H}_2\text{C}-\text{C}-\text{CH}_2-\text{O}-\text{NO}_2 \\   \\ \text{CH}_2 \\   \\ \text{O} \\   \\ \text{NO}_2 \end{array}$	316,2	141,3	202-05	
4	Nitrocellu- lose	NC	Bomullskrut Cellulosanitrat		~ 10 <sup>5</sup>	-	185-190	
5	Hexogen	RDX	Cyclonit, T4, 1,3,5-Trinitro- s-triazin, 1,3,5-trinitro- cyklohexan		222,1	202	230	Decomposition close above the melting point

contin. Table 1

No	Explosive	Abbrev.	Synonym	Formula	Mol. weight	Melt. point	Deflagr. temp. °C	Remark
6	Octogen	HMX	Homocyclonit, 1,3,5,7-Tetra- nitro-1,3,5,7- tetraazacyklo- oktan		296,2	282	290	Modifications: (α orthoromb) (β monocline) (γ monocline) (δ hexagonal)
7	Tetryl	TETR	CE, 2,4,6,N-Tetra- nitro-N-metyl- anilin		287,2	131,5	185-195	Melts with decomposition
8	Ammonium picrate	AM-PIKR	Explosive D, Ammonium- 2,4,6-trinitro- rofenolat		246,1	265-71	320	Melts with decomposition
9	Hexanitrostilbene	HNS	2,2',4,4',6,6'- Hexanitrostil- ben		450	316		Decomposition at 280 °C: 5h 7 % 7h 26 % 8h 48 %

contin. Table 1














No	Explosive	Abbrev	Synonym	Formula	Mol.- weight	Melt.- point	Deflagg. temp. °C	Remark
10	Trinitrobenzene	TNB 1,3,5-TNB	Benzit, 1,3,5-Trinitrobenzen		213,1	123,2		
11	Trotyl, 2,4,6-Trinitrotoluene	TNT 2,4,6-TNT	2,4,6-Trinitrotoluen Tri, Trotyl, Trolit, sym-TNT, α-TNT		227,1	80,8	300	Decomposition begins at 180 °C



Table 2. List of nitrotoluenes\*. Formulas and physical data.

∞

No	Compound	Abbrev	Synonym	Formula	Mol.-weight	Melt.-point
12	2,4,5-Trinitrotoluene	2,4,5-TNT	γ-TNT		227,1	104
13	2,3,6-	2,3,6-TNT	η-TNT			111
14	2,3,5-	2,3,5-TNT	ε-TNT			97,2
15	2,3,4-	2,3,4-TNT	β-TNT			112
16	3,4,5-	3,4,5-TNT	δ-TNT			137,5
17	2,6-Dinitrotoluene	2,6-DNT			182,1	66
18	2,5-	2,5-DNT				52,5
19	2,4-	2,4-DNT				70,1
20	2,3-	2,3-DNT				61
21	3,5-	3,5-DNT				93
22	3,4-	3,4-DNT				59,8

\*2,4,6-Trinitrotoluene, see no 11, table 1.

Contin. Table 2




No	Compound	Abbrev	Synonym	Formula	Mol. weight	Melt. point
23	2-Mononitrotoluene	2-MNT	o-nitro-toluen		137, 1	{ -10, 6(α) - 4, 1(β)
24	"	3-MNT	m-		24	15, 5
25	"	4-MNT	p-		25	51, 3

Table 3. Solubilities of explosives (g/100 g at 25 °C unless specified)

Solvent Explosive	Acetone	Aceto- nitrile	Benzene	Butyro- lactone	Dimethyl- formamide	Dimethyl- sulph- oxide	Ethanol	Ether	Ethyl- acetate	Chloro- form	Carbon- disul- phide	Carbon- tetra- chloride	Toluene	Water
1. EGDN	∞		∞				v s	∞	∞	∞	sl s	ca 2 vol%	∞	0,56
2. NG	∞		∞				54 (20°C)	∞	∞	∞	1,25 vol%	2 vol%	v s	0,15 (20°C)
3. PETN	20,3 (20°C)		0,3 (20°C)				0,20 (20°C)	0,25 (20°C)	6,3 (19°C)	0,06 (19°C)			0,23 (20°C)	sl s
4. NC <sup>1</sup>	v s													
5. RDX	8,2	5,5	0,05 (20°C)	14	37	41	0,11 (20°C)	0,055 (20°C)	1,5 (20°C)	0,008 (20°C)	sl s	0,005 (50°C)	0,02 (20°C)	0,006
6. HMX	2,8	2,0		12	2,3 <sup>2</sup> (20°C)	57								sl s
7. TETR	68 vol%		3,5 (20°C)		114		0,65	0,46	12,2 (18°C)	0,68	0,024	0,031	3,0 (19,5°C)	0,008
8. AM-PIKR					90	75	0,62 (24°C)	0,004 (24°C)						1,1 (20°C)
9. HHS	< 0,1			0,4	1,5	1,4	s							sl s
10. TNB(17°C)	59		6,2				2,1	1,7	30	6,2	0,24	0,24	12	0,03 (15°C)
11. TNT	132		88		142	128	1,5	3,8		25	0,63	0,82	67	0,02

1) Ketones, esters and ether + ethanol are common solvents for nitrated cellulose. The solubility depends on degree of nitration.

2) Solvate crystals precipitate.

## CHROMATOGRAPHY

### High pressure liquid chromatography (HPLC)

High pressure liquid chromatography is used in qualitative as well as in quantitative analysis<sup>2</sup>. Substances of low vapour pressure or of thermal instability are often more easy to analyse by means of high pressure liquid chromatography than by gas chromatography.

The instrumentation used was a Waters Assoc liquid chromatograph, model ALC-GPC-204, equipped with a 254 UV detector, model 440. A  $\mu$ -Porasil column was used (1/4" x 1", 10  $\mu$ , 400 m<sup>2</sup>/g, Waters Assoc). Eluents were of Uvasol (Merck's) purity. The sensitivity of the detector is approximately of the same order for the aromatic nitro compounds and the nitramines but some ten powers lower in the case of the nitric esters<sup>1</sup>.

Retention volumes of the various explosives are given in table 4. A solvent peak emanating from the sample injected may appear in a HPLC chromatogram. In case the peak is of considerable height it has been indicated.



Reversed phase chromatography (Lichrosorb C2 column and water-ethanol as eluent) resulted in chromatograms with less resolved peaks.

Table 4. Retention volumes (ml) of explosives on a  $\mu$ -Porasil column.  $V_M = 3,0$  ml (t = tailing).

Eluent Substance	1	2	3	4	5
11. TNT	3,1	3,1	3,3	5,6	5,0
1. EGDN	3,2	3,3	3,8	4,8	6,1
2. NG	3,2	3,4	4,0	5,2	7,0
10. TNB	3,1	3,3	4,0	7,3	7,5
3. PETN	3,2	3,2	7,0	15,6	8,6
7. TETR	3,3	3,6	7,0	15,9	31,0
9. HNS	3,2	3,2	3,6	> 40	> 40
5. RDX	4,8	11,4	> 40	> 40	
6. HMX	8,6	> 40	> 40		
CHCl <sub>3</sub>			3,0	3,3	3,2
Acetone	3,7	3,8	3,9	t	t

Eluents:

1. Chloroform-acetonitrile 90:10
2. Chloroform
3. c-Hexane-chloroform 50:50
4. c-Hexane-dichloromethane 70:30
5. c-Hexane-chloroform 80:20



Column:  
μ-Porasil

Eluent:  
c-Hexane-dichloromethane, 70:30

Flow rate:  
2,0 ml/min

Sample volume:  
10 μl

Sample conc:  
Totally 2 % in ethanol  
EGDN:NG 1:1


Column:  
μ-Porasil

Eluent:  
Chloroform-acetonitrile, 90:10


Flow rate:  
1,0 ml/min

Sample volume:  
3 μl


Sample conc:  
Totally 0,6 % in acetone  
TNT:RDX:HMX 1:1:1




Column:  
 $\mu$ -Porasil  
 Eluent:  
 Chloroform  
 Flow rate:  
 1,0 ml/min  
 Sample volume:  
 10  $\mu$ l  
 Sample conc:  
 Totally 0,07 % in chloroform  
 TNT:TETR 3:4



Column:  
 $\mu$ -Porasil  
 Eluent:  
 c-Hexane-chloroform,  
 80:20  
 Flow rate:  
 2,0 ml/min  
 Sample volume:  
 20  $\mu$ l  
 Sample conc:  
 Totally 0,03 % in  
 chloroform  
 TNT:PETN 1:300



Column:  
 $\mu$ -Porasil  
 Eluent:  
 c-Hexan-chloroform,  
 80:20  
 Flow rate:  
 1,0 ml/min  
 Sample volume:  
 20  $\mu$ l  
 Sample conc:  
 Totally 0,01 % in  
 chloroform  
 TNT:TNB 1:1




Column:  
 $\mu$ -Porasil

Eluent:  
c-Hexan-dichloromethane, 87,5:12,5


Flow rate:  
2,0 ml/min

Sample volume:  
20  $\mu$ l

Sample conc:  
Totally 0,02 %



Column:  
 $\mu$ -Porasil  
 Eluent:  
 c-Hexan-chloroform, 95:5  
 Flow rate:  
 2,5 ml/min  
 Sample volume:  
 10  $\mu$ l  
 Sample conc:  
 Totally 0,1 %



Column:  
 $\mu$ -Porasil  
 Eluent:  
 c-Hexane  
 Flow rate:  
 1,0 ml/min  
 Sample volume:  
 10  $\mu$ l  
 Sample conc:  
 Totally 0,05 %



Gas chromatography (GC)

Gas chromatography is used in qualitative as well as in quantitative analysis and has been applied in analyses of hexoto1 samples<sup>8</sup>.

The instrument used was a Varian gas chromatograph model 3700, equipped with a flame ionization detector. Columns of Pyrex glass were used, having OD 1/4", ID 1.5 mm, and the lengths of 0.75 m and 3 m. Data on column length, stationary phase, flow rates of gases, temperatures and sample size are given below each chromatogram. In all cases a sample volume of 1  $\mu$ l was injected.


Because of their low vapour pressures the substances Nos 4(NC), 6(HMX), 8(AM-PIKR) and 9(HNS) were not chromatographed.

1 EGDN

3 PETN

17

2 NG



Column: 5% DC 550 on Chromosorb W-HP, 0,75 m x 1,5 mm

Temp:

Column: Progr Start: 75°C during 4 min, 10°/min  
Final temp 200 °C

Inj: 130°C

Det: 220°C

Flow rate:

N<sub>2</sub> 20 ml/min


H<sub>2</sub> 20 ml/min

Air 250 ml/min

Sample conc: Totally 2 % in acetone

EGDN:NG:PETN = 1:2:3

18    11    TNT  
       5    RDX




Column: 5 % DC 550 on  
 Chromosorb W-HP  
 0,75 m x 1,5 mm

Temp:  
 Column: 160°C isotherm  
 Inj: 170°C  
 Det: 190°C

Flow rate:  
 N<sub>2</sub> 20 ml/min  
 H<sub>2</sub> 20 ml/min  
 Air 250 ml/min

Sample conc:  
 Totally 0,35 % in acetone  
 TNT : RDX = 2 : 7

11    TNT  
       7    TETR



Column: 5 % DC 550 on  
 Chromosorb W-HP  
 0,75 m x 1,5 mm

Temp:  
 Column: 190°C isotherm  
 Inj: 190°C  
 Det: 200°C

Flow rate:  
 N<sub>2</sub> 20 ml/min  
 H<sub>2</sub> 20 ml/min  
 Air 250 ml/min

Sample conc:  
 Totally 0,5 % in acetone  
 TNT : TETRYL = 1 : 2


10 TNB

16 - 22 DNT

11 - 15 TNT

23 - 25 MNT

19



Column: 3 % OV 225 on Chromosorb W-HP, 3 m x 1,5 mm

Temp:

Column: Progr Start: 130°C during 4 min, 10°/min  
Final temp 220 °C

Inj: 150°C

Det: 230°C

Flow rate:

N<sub>2</sub> 20 ml/min

H<sub>2</sub> 20 ml/min

Air 250 ml/min

Sample conc: Totally 0,5 % in acetone

MNT:DNT:TNT:TNB = 3:12:6:1

## SPECTROSCOPY

The purity of the samples was verified by gas chromatography.

### Mass spectrometry (MS)

The mass spectra were recorded with a LKB Model 2091 single focus <sup>9</sup> spectrometer. The substances were analyzed by EI (electron impact) at electron voltages of 70 and 20 eV and by CI (chemical ionization) with methane (1 Torr) as ionization media.

The majority of the samples were introduced into the mass spectrometer via a gas chromatograph. The columns used were identical to those specified in the chapter of gas chromatography. The substances Nos. 6(HMX) and 7(TETR) were introduced through the direct inlet, those of low vapour pressure Nos. 4(NC), 8(AM-PIKR), and 9(HNS) were not examined.

Data on mode of ionization, electron voltage, ion source temperature and characteristic fragments are given together with the spectrum.

### Infrared spectroscopy (IR)

The IR spectra were recorded with a Perkin Elmer grating spectrometer model 377. Time constant: Auto; Scan mode: 13 min; Slit: N.

The majority of the samples were prepared in the form of a KBr disc composed of 0.5 - 1.8 mg substance in 300 mg KBr (diameter 13 mm, thickness 1 mm, applied pressure 1.500 kg/cm<sup>2</sup>). The substances Nos. 1, 2, and 23 - 25 were dissolved in CCl<sub>4</sub> (0.1 - 0.2 M) and analyzed in a 0.12 mm NaCl-cell. Substance No. 4(NC) was analyzed as a film prepared from a solution of the sample in acetone.

Significant band assignments and sample concentration are given at each spectrum.

Table 5. IR-frequencies of (x)-NO<sub>2</sub> associated with nitric esters, nitramines and nitro compounds.

Substance group	(x)-NO <sub>2</sub> (cm <sup>-1</sup> )	X
Nitric esters Nos. 1 - 4	1660 - 1640; 1285 - 1270	O
Nitramines Nos. 5 - 6	1590 - 1530; 1310 - 1270	N
Nitro compounds Nos. 7 - 25	1560 - 1520; 1370 - 1340	C

## Nuclear magnetic resonance (NMR)

This method is used both in quantitative and qualitative analysis and has been applied in analyses of hexitol samples<sup>13</sup>.

### <sup>1</sup>H NMR

Proton spectra<sup>14</sup> were obtained at 60 MHz on a Varian NV-14 spectrometer (CW). The internal reference was tetramethylsilane (TMS). The spectra were recorded with a sweep time of 500 seconds. Sample concentrations were 2 - 3 % in acetone-d<sub>6</sub>\*, except 9(HNS) and 4(NC) which were dissolved in dimethylsulfoxide-d<sub>6</sub>\*\* and 1(EGDN) dissolved in CCl<sub>4</sub>. To 2(NG) acetone was added.

The chemical shifts of the different protons are indicated by letter symbols in formula and spectrum. Certain parts of the spectra have been expanded.

### <sup>13</sup>C NMR

Carbon-13 spectra<sup>15</sup> were obtained at 20 MHz with a Varian CFT 20 spectrometer system, equipped with a 10 mm probe. The pulse interval varied between 1 to 3 seconds. The internal reference was tetramethylsilane (TMS) and the spectra obtained are proton decoupled.


Sample concentrations were 3 - 8 % w/v in acetone-d<sub>6</sub>, except 2(NG) and 9(HNS) which were 5 % in CCl<sub>4</sub> and 1 % in dimethylsulfoxide-d<sub>6</sub>, respectively. To shorten the relaxation times, Cr(acac)<sub>3</sub> was added to a concentration of  $3.6 \times 10^{-2}$  M, except for 1 (EGDN) and 4(NC).

The line assignments are indicated by a letter at the corresponding carbon atom in the structure. The large solvent peaks from acetone-d<sub>6</sub> and dimethylsulfoxide-d<sub>6</sub> around 30 and 40 ppm, respectively, are not presented.


\* Signal at 2.05 ppm

\*\* Signal at 2.50 ppm


## EGDN




Inlet: GC  
Column: DC 550  
Ion source: 150 °C



Inlet: GC  
Column: DC 550  
Ion source: 150 °C




Inlet: GC  
Column: DC 550  
Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ )

2965, 2910, 2920	C-H aliph	1460, 1428	C-H aliph	837	O-( $\text{NO}_2$ )
1655	(O)- $\text{NO}_2$ asym	1270	(O)- $\text{NO}_2$ sym		




$^1\text{H}$  60 MHz




$^{13}\text{C}$  20 MHz









Group frequencies ( $\text{cm}^{-1}$ )



2961, 2923 C-H aliph

1424 C-H aliph


830 O-( $\text{NO}_2$ )

1660 (O)- $\text{NO}_2$  asym


1272 (O)- $\text{NO}_2$  sym




26  
PETN




Inlet: GC  
Column: DC 550  
Ion source: 150 °C



Inlet: GC  
Column: DC 550  
Ion source: 150 °C




Inlet: GC  
Column: DC 550  
Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ )

3280	Over tone from 1640	1660, 1640	(O)-NO <sub>2</sub> asym	1286, 1278	(O)-NO <sub>2</sub> sym	702	(O)-NO <sub>2</sub>
3024, 2987	C-H aliph	1470, 1397	C-H aliph	870, 850	O-(NO <sub>2</sub> )		

<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz



28

NC


See page 20.




Group frequencies ( $\text{cm}^{-1}$ )

3400 OH free	3020, 2940 C-H aliph	1440 C-H aliph	1082 C-O-C	765, 705, 550 (O)-NO <sub>2</sub>
3350 OH bound	1670, 1650 (O)-NO <sub>2</sub> asym	1295 (O)-NO <sub>2</sub> sym	855 O-(NO <sub>2</sub> )	


<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




RDX




Inlet: GC  
 Column: DC 550  
 Ion source: 150 °C



Inlet: GC  
 Column: DC 550  
 Ion source: 150 °C




Inlet: GC  
 Column: DC 550  
 Ion source: 150 °C




Group Frequencies ( $\text{cm}^{-1}$ )

3079, 3069, 3002, 2951	C-H aliph	1459, 1389	C-H aliph	1040, 923	C-H
1590, 1570, 1531	(N)-NO <sub>2</sub> asym	1310, 1269	(N)-NO <sub>2</sub> sym	702, 602	(N)-NO <sub>2</sub>






## HMX



Inlet: Directly

Temp: 105 °C


Ion source: 150 °C



Inlet: Directly

Temp: 105 °C


Ion source: 150 °C



Inlet: Directly

Temp: 105 °C

Ion source: 120 °C



Group frequencies ( $\text{cm}^{-1}$ ):

3042, 3032, 2986 C-H aliph


1462, 1395 C-H aliph

1147, 968 Ringabsorption


1560, 1540, 1542 (N)-NO<sub>2</sub> asym

1298, 1280 (N)-NO<sub>2</sub> sym







Inlet: Directly  
Temp: 85 °C  
Ion source: 150 °C



Inlet: Directly  
Temp: 85 °C  
Ion source: 150 °C



Inlet: Directly  
Temp: 85 °C  
Ion source: 125 °C




Group frequencies (cm<sup>-1</sup>):

3100, 3084, 3022	C-H arom	1611	phenyl	1354, 1341	(C)-NO <sub>2</sub> sym
2958	C-H aliph	1552, 1545, 1539	(C)-NO <sub>2</sub> asym och (N)-NO <sub>2</sub>	1287	(N)-NO <sub>2</sub> sym
				1080, 922	Isol C-H arom

<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz



36



AM-PIKR

See page 20



Group frequencies ( $\text{cm}^{-1}$ ):


3300, 3250	$\text{NH}_4^+$	3020	$\text{NH}_4^+$	1570, 1545	(C)- $\text{NO}_2$ asym	1333	(C)- $\text{NO}_2$ sym
3100	C-H arom	1624, 1604	phenyl	1422	$\text{NH}_4^+$	1080, 910	isol C-H arom



38



HNS

See page 20




Group frequencies (cm<sup>-1</sup>):

3105, 3083, 3074	C-H arom	1543, 1537 (C)-NO <sub>2</sub> asym	960	Transalken
1618, 1611	phenyl	1353, 1349 (C)-NO <sub>2</sub> sym	1087, 922	C-H arom







## TNB; 1,3,5-TNB




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C





Group frequencies (cm<sup>-1</sup>):


3115, 3052 C-H arom  
1622, 1440 phenyl

1542 (C)-NO<sub>2</sub> asym  
1347 (C)-NO<sub>2</sub> sym


1080, 922 Isol C-H arom




## TNT: 2,4,6-TNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C





Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):


3104, 3060, 3021	C-H arom	1622, 1604	phenyl	1358	(C)-NO <sub>2</sub> sym
2960	C-H aliph	1544, 1536	(C)-NO <sub>2</sub> asym	1090, 910	isol C-H arom




## 2,4,5-TNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ ):

3120, 3060 C-H arom	1563, 1547, 1537 (C)-NO <sub>2</sub> asym	1040, 918 Isol C-H arom
1612, 1604 phenyl	1350 (C)-NO <sub>2</sub> sym	

<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




a. 20.17  
b. 123.27  
c. 130.73  
d. 140.83  
e. 142.48  
f. 145.02  
g. 151.18


## 2,3,6 - TNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):


3108 C-H arom 1367, 1345 (C)-NO<sub>2</sub> sym

1560, 1511, 1532 (C)-NO<sub>2</sub> asym

<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz







## 2,3,5-TNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):

3108, 3064 C-H arom      1563, 1545 (C)-NO<sub>2</sub> asym      1080, 919 Isol C-H arom  
1630, 1599 phenyl      1352, 1348 (C)-NO<sub>2</sub> sym

<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




a. 17.45  
b. 120.08  
c. 133.26  
d. 136.23  
e. 141.73  
f. 147.45  
g. 148.62


## 2,3,4-TNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group Frequencies (cm<sup>-1</sup>):

3097, 3060 C-H arom 1560, 1553, 1548 (C)-NO<sub>2</sub> asym  
1608, 1597 phenyl 1357, 1348 (C)-NO<sub>2</sub> sym


<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




## 3,4,5 - TNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):

3082, 3060 C-H arom	1370, 1350, 1342 (C)-NO <sub>2</sub> sym
1563, 1542 (C)-NO <sub>2</sub> asym	929 Isol C-H arom

<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




a. 21.33  
b. 137.12  
c. 136.38  
d. 142.06  
e. 146.58


## 2,6-DNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C





Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ ):

3108, 3100, 3038	C-H arom	1614	phenyl	1365, 1352	(C)-NO <sub>2</sub> sym
3000, 2957	C-H aliph	1542, 1530, 1527	(C)-NO <sub>2</sub> asym		







## 2,5 - DNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):

3124, 3100, 3044	C-H arom	1547, 1532, 1523	(C)-NO <sub>2</sub> asym	903	Isol C-H arom
2992	C-H aliph	1350	(C)-NO <sub>2</sub> sym		


<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




58  
2,4-DNT




Inlet: GC  
Column: OV 225  
Ion source: 150 °C



Inlet: GC  
Column: OV 225  
Ion source: 150 °C




Inlet: GC  
Column: OV 225  
Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):

3105, 3085, 3060	C-H arom	1540, 1523, 1509	(C)-NO <sub>2</sub> asym	1070, 913	Isol C-H arom
1608	phenyl	1348	(C)-NO <sub>2</sub> sym		


<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




60  
2,3-DNT




Inlet: GC  
Column: OV 225  
Ion source: 150 °C



Inlet: GC  
Column: OV 225  
Ion source: 150 °C




Inlet: GC  
Column: OV 225  
Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ ):

3105 C-H arom 1547, 1535 (C)-NO<sub>2</sub> asym  
1610 phenyl 1370, 1350 (C)-NO<sub>2</sub> sym


<sup>1</sup>H 60 MHz




<sup>13</sup>C 20 MHz




## 3,5-DNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C





Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):

3120, 3088, 3074	C-H arom	1542, 1531, 1523	(C)-NO <sub>2</sub> asym	1078, 918, 902	isol C-H arom
1590	phenyl	1353, 1341	(C)-NO <sub>2</sub> sym		







## 3,4 - DNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C





Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ ):


3118, 3078, 3061	C-H arom	1544, 1527	(C)-NO <sub>2</sub> asym	899	isol C-H arom
1608, 1601	phenyl	1364, 1350	(C)-NO <sub>2</sub> sym		




66  
2-MNT




Inlet: GC  
Column: OV 225  
Ion source: 150 °C



Inlet: GC  
Column: OV 225  
Ion source: 150 °C





Inlet: GC  
Column: OV 225  
Ion source: 150 °C




Group frequencies (cm<sup>-1</sup>):


3078	C-H arom	1615	phenyl	1528	(C)-NO <sub>2</sub> asym
2984, 2940	C-H aliph	1580	C-H arom	1351	(C)-NO <sub>2</sub> sym




## 3-MNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C





Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ ):


3100, 3078, 3039	C-H arom	1590	phenyl	1354	(C)-NO <sub>2</sub> sym
2932	C-H aliph	1532	(C)-NO <sub>2</sub> asym		




## 4-MNT




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C



Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Inlet: GC  
 Column: OV 225  
 Ion source: 150 °C




Group frequencies ( $\text{cm}^{-1}$ ):

3085, 3053	C-H arom	1608, 1601	phenyl	1523	(C)- $\text{NO}_2$ asym
2992, 2956, 2936	C-H aliph	1570	C-H arom	1350	(C)- $\text{NO}_2$ sym

$^1\text{H}$  60 MHz



$^{13}\text{C}$  20 MHz





## COLOUR REACTIONS

### Spot tests<sup>16</sup>

The reactions were performed using a sample size of 10 µg dissolved in 50 µl acetone (0.02 %). Two drops of each reactant were added. Ethanol can be chosen instead of acetone if the sample is soluble in this solvent.


### Methods of detection

1. Sodium hydroxide  
1 M NaOH  
Colour appears instantly
2. Sodium hydroxide + Griess' reagent  
A. 1 M NaOH  
B. Griess' reagent:  
1 vol 1 % sulphanilic acid in 30 % HAc + 1 vol 0.1 %  
N-(1-naphthyl)-ethylenediaminodihydrochloride in water  
10 min between adding A and B, colour appears in 2 min.
3. Zinc + Griess' reagent  
A. Zinc dust  
B. Griess' reagent as above  
Colour appears in 2 min.
4. Zinc + hydrochloric acid + p-dimethylaminobenzaldehyde (DMAB)  
A. Zinc dust  
B. 3 M HCl  
C. Saturated solution of DMAB in benzene.  
Add A + B, boil for 2 min and centrifugate. Transfer the solution to a filter paper and dry. Add C. Colour appears instantly (heating may be necessary).
5. Diphenylamine (DPA) in sulphuric acid  
0.1 % DPA in conc. H<sub>2</sub>SO<sub>4</sub>  
Colour appears instantly.

Reactions:


Method 1: Polynitroaromates are transformed to coloured complexes in alkaline media (Meisenheimer complex).

Method 2: Nitrite, liberated by alkali, and Griess' reagent react giving an azo dyestuff as product (within 2 min).



Method 3: Nitrate, liberated by acid, is reduced to nitrite by Zn + acid. The colour reaction is the same as in Method 2.

Method 4: Nitrosubstituents of aromates is reduced to primary amines by Zn + acid. These react with DMAB giving coloured Schiff's bases as products.



Method 5: Nitrate or nitrite, liberated by acid, oxidizes DPA to a blue dyestuff.

Table 6. Colours with spot test reagents.

Detection method Substance	1	2	3	4	5
1. EGDN	-	Or +	Or ++	-	B1 ++
2. NG	-	Or +++	-"-	-	B1 +++
3. PETN	-	Or + <sup>1</sup>	-"-	-	B1 +++ <sup>1</sup>
4. NC	-	Or +++	Rö + <sup>1</sup>	-	B1 ++ <sup>1</sup>
5. RDX	-	Or +++	Or ++	-	B1 +++ <sup>1</sup>
6. HMX	-	Or ++	Or +	-	- <sup>1</sup>
7. TETR	RöBr +++	RöOr +++ <sup>1</sup>	- <sup>1</sup>	Or +	- <sup>1</sup>
8. AM-PIKR	Gu ++	GuGr ++	-	Gr +	-
9. HNS	RöBr +++	BrGu +	-	Or +	-
10. 1,3,5-TNB	Br +++	Gu ++ <sup>1</sup>	-	Gu +++	-
11. 2,4,6-TNT	RöBr +++	GuBr ++ <sup>1</sup>	-	GuOr ++	-
12. 2,4,5-TNT	RöVi +++	Or ++	-	RöVi ++	-
13. 2,3,6-TNT	BrGu ++	-"-	-	Gu +	-
14. 2,3,5-TNT	Rö ++	-"-	-	Or ++	-
15. 2,3,4-TNT	RöBr +++	-"-	-	Gu +	-
16. 3,4,5-TNT	Vi +++	-"-	-	Gu ++	-
17. 2,6-DNT	Rö + <sup>1</sup>	- <sup>1</sup>	-	Gu +++	-
18. 2,5-DNT	- <sup>1</sup>	- <sup>1</sup>	-	Or +++	-
19. 2,4-DNT	B1 ++ <sup>1 2</sup>	Rö ++ <sup>1</sup>	-	OrGu +++	-
20. 2,3-DNT	- <sup>1</sup>	- <sup>1</sup>	-	Gu +	-
21. 3,5-DNT	Rö + <sup>1</sup>	Rö ++ <sup>1</sup>	-	GuOr +++	-
22. 3,4-DNT	- <sup>1</sup>	- <sup>1</sup>	-	Gu ++	-
23-25. MNT-er	- <sup>1</sup>	- <sup>1</sup>	-	Gu +++	-
Nitrite	-	Rö +++	Rö +++	-	B1 +++
Nitrate	-	-	Rö +++	-	B1 +++

<sup>1</sup> Valid for 100 µg substance<sup>2</sup> Colour fades rapidly

## Abbreviations:

B1	blue	Rö	red	-	no colour
Gr	green	Vi	violet	+	faint
Gu	yellow	Br	brown	++	medium
Or	orange			+++	strong

### Detection on silica gel layers

The data of reagents and colours given below are cited from a report of Bilson<sup>17</sup>. These reagents have similar underlying reaction mechanisms as the spot tests. The modifications of composition of the reagents, conditioned by the presence of silica gel, will in some cases cause a change in colour shade.

Spotted amounts of sample have been 10 µg. Colours were developed by spray techniques.

### Methods of detection

1. No reagents
  - a) Daylight
  - b) UV-light, 254 nm
2. Alkaline methanol  
10 % KOH in methanol
3. Alkaline methanol + Griess' reagent
  - A. Alkaline methanol, see above
  - B. Griess' reagent, see spot test reagent 2 B.20 min at 105 °C between A and B.
4. Diphenylamine (DPA)  
1.5 % DPA in methanol  
UV-light, 254 nm, 15 min.
5. Titaniumtrichloride + p-dimethylaminobenzaldehyde (DMAB)
  - A. 12.5 %  $\text{TiCl}_3$  in 15 % HCl
  - B. 0.25 % DMAB<sup>3</sup> in 30 % HAc
  - a) Air drying between A and B
  - b) + heating at 105 °C during 15 min.

Table 7. Colours on silica gel layers (cit. Bilson<sup>17</sup>).

Detection method Substance	1a	1b	2	3	4	5a	5b
2. NG	-	-	-	Rö ++	GBr ++	-	-
3. PETN	-	-	-	""	""	-	-
4. NC	-	Pu (+)	-	""	GuGr +	-	-
5. RDX	-	Pu ++	-	""	31G ++	Gu ++	Gu ++
6. HMX	-	""	-	""	Gu +	-	-
7. TETR	Gu ++	BrPu ++	RöBr ++	GuBr ++	Br ++	GuOr ++	Br ++
8. AM-PIKR	""	""	Gu ++	Gu ++	Gu ++	Gu ++	OrBr ++
10. TNB	-	Pu ++	RöBr ++	RöBr +	PuBr ++	GuOr ++	GuBr ++
11. TNT	-	""	""	RöBr ++	Br +	Gu ++	""
19. 2,4-DNT	-	""	-	RöBr +	Br +	""	Gu ++
23. 2-MNT	-	""	-	-	-	PuBr ++	PuBr ++

## Abbreviations:

B1	blue	Gu	yellow	-	no colour
Br	brown	Or	orange	(+)	very pale
G	grey	Pu	purple	+	pale
Gr	green	Rö	red	++	distinct

#### THIN LAYER CHROMATOGRAPHY

Thin layer chromatography is used in quantitative and qualitative analysis<sup>2</sup>.

The substances were chromatographed on 0.25 mm silica gel layer containing fluorescence indicator (Silica Gel F-254, Merck No. 5715). Spotted amount was 2  $\mu$ l of a 0.1 % solution (2  $\mu$ g) with the exception of the nitric esters where sample amount was 10 - 20 times bigger.

The eluents were chosen according to Bilson<sup>17</sup>. The experimental conditions regarding activity of layers and saturation of chamber, however, were different. Found and cited  $R_f$ -values in table 8 therefore correspond to those maximum and minimum values one can expect for silica layers. The travel distance of solvents were 10 cm.

One example of a two dimensional chromatogram is given in figure 1.




Figure 1. Twodimensional chromatogram of explosives (no 1 - 3, 5 - 11). Layer: Silica Gel F-254. Travel distance: 10 cm.

Direction 1: Light petroleum (60 - 80 °C)/Diisopropylether 25:75

Direction 2: Benzene/Acetone 70:30

Table 3.  $R_f$ -values on silica gel layers.

A. Layers not activated. Chamber not lined with filterpaper.

Substances	Eluents										Two dimen	
	1	2	3	4	5	6	7	8	9	10	1	2
3 AM-PIKR	0	0	0	0	0,11	0,14	0,38	0,36	0,74	0,97	0	0,12
6 HMX	0	0	0,01	0,02	0,36	0,25	0,56t	0,66	0,85	0,76u	0	0,43
5 RDX	0	0,01	0,03	0,10	0,48	0,72	0,64	0,71	0,80	0,81	0,01	0,53
9 HNS	0	0	0	0,22	0,95	0,80	0,96	0,97	0,99	0,95	0,02	0,93
7 TETR	0	0,03	0,18	0,43	0,88	0,82	0,89	0,94	0,98	0,92	0,17	0,87
2 NG	0,01	0,12	0,43	0,66	0,87	0,91	0,88	0,93	0,98	0,92	0,42	0,87
3 PETN	0	0,13	0,66	0,66	0,93	0,91	0,93	0,94	0,99	0,94	0,55	0,92
1 EGDN	0,03	0,23	0,42	0,73	0,87	0,92	0,90	0,92	0,97	0,93	0,58	0,85
10 TNB	0	0,26	0,58	0,63	0,94	0,87	0,94	0,96	0,99	0,93	0,70	0,92
11 TNT	0	0,45	0,71	0,73	0,96	0,90	0,95	0,96	0,99	0,94	0,87	0,92

Eluents in twodimensional chromatogram: Direction 1. Light petroleum (60 - 80 °C)/Diisopropylether 25:75  
 Direction 2. Benzene/Acetone 70:30



atin. Table 3.

B. Layers activated 1 hour at 105 °C, chamber lined with filter paper (cit Bilson<sup>17</sup>).

Eluent Substance	1	2	3	4	5	6	7	8	9	10
8 AM-PIKR	0	0	0	0	0,13	0,17	0,25	0,29	0,53	0,80
6 HMX	0	0	0	0	0,31	u	u	0,46	0,59	u
5 RDX	0	0	0,02	0,07	0,42	0,52	0,48	0,53	0,57	0,64
9 HNS	-	-	-	-	-	-	-	-	-	-
7 TETR	0	0,03	0,13	0,33	0,71	0,70	0,72	0,75	0,77	0,77
2 NG	0	0,06	0,31	0,47	0,70	0,82	0,70	0,71	0,76	0,78
3 PETN	0	0,08	0,51	0,49	0,75	0,84	0,76	0,79	0,76	0,77
1 EGDN	-	-	-	-	-	-	-	-	-	-
10 TNB	0	0,17	0,42	0,47	0,76	0,81	0,77	0,81	0,78	0,76
11 TNT	0	0,28	0,52	0,54	0,76	0,83	0,79	0,79	0,78	0,76

t = tailing u = 1000

Eluents (mixtures in v/v)

1. Light petroleum, 60-80 °C
2. Light petroleum, 60-80 °C/Diisopropylether 70:30
3. Light petroleum, 60-80 °C/n-Butyl acetate 80:20
4. Benzene
5. Benzene/Acetone 75:25
6. Dichloromethane/Acetic acid 90:10
7. Benzene/Methyl alcohol/Acetic acid 90:16:8
8. Benzene/Acetone/Methyl alcohol 70:20:10
9. Dioxan/benzene/Acetic acid 75:27:12
10. Methyl alcohol

Purity:

- Acetone pa  
Benzene pa  
n-Butyl acetate, 98-100 °C, pure  
Diisopropylether, 66-69 °C, ca 97  
Dichloromethane, 98-100  
Dioxan pa  
Methyl alcohol pa  
Light petroleum, 60-80 °C, pa  
Acetic acid, 99-100 %, pa
- Merck no 14  
Merck no 1783  
Kebo no 25477  
Kebo no 23723  
Kebo no 13090  
Merck no 9671  
Merck no 6009  
Merck no 1774  
Riedel-de Haen  
no 32209

## REFERENCES

1. Alm, A: Detection and analysis of explosives. FOA Report C 20133-D1. Aug. 1976. (In Swedish)
2. Alm, A: Explosives, gunpowders and related compounds. I. Literature list on liquid chromatographic separations. II. Literature list on gas chromatographic separations. FOA Reports C 20101-D1 and C 20102-D2. March 1976 (In Swedish)
3. Hogan, V D, Richter, T A E: A New Convenient Tool for Identifying Composite Explosives: Proton Magnetic Resonance Fingerprinting. Technical Report 4790, Picatinny Arsenal, Dover, New Jersey, June 1975.
4. Fedoroff, B T (Ed.): Encyclopedia of explosives and related items. (Picatinny Arsenal, USA). (Vol. 7 published 1975).
5. Meyer, R: Explosivstoffe. 4 Ed. Weinheim 1975.
6. Analytical methods for powders and explosives. Ed. AB Bofors Nobelkrut. 1960.
7. Schroeder, W A et al.: Ultraviolet and visible absorption spectra in ethyl alcohol. Anal. Chem. 23 (1951) 1740-47.
8. Frölén-Lindgren, I: Quantitative analysis of hexotol with gas chromatography. FOA Report C 20178-D1. June 1977. (In Swedish)
9. Volk, F, Schubert, H: Massenspektrometrische Untersuchungen von Explosivstoffen. Explosivstoffe 16 (1968) 2-10.
10. Gillis, R G, Lacey, M J, Shannon, J S: Chemical ionization mass spectra of explosives. Org. Mass Spectrom. 9 (1974) 359-64.
11. Pristera, F, Fredericks, W: Compilation of infrared spectra of ingredients of propellants and explosives. Picatinny Arsenal Technical Memorandum 1887. Sept. 1969 (AD 859 846), see also Pristera et al. Anal. Chem. 32 (1960) 495-508, 25 (1953) 844-56.
12. Bellamy, L J: The infra-red spectra of complex molecules. 2 Ed. London 1964.
13. Alm, A, Östman, B: Quantitative analysis of hexotol with nuclear magnetic resonance. FOA Report C 20015-D1. Oct. 1974. (In Swedish)
14. Varian Associates High Resolution NMR Spectra Catalog. Combined Vols. 1 and 2. N S Bhacca, D P Hollis, L F Johnson, E A Pier, and J N Shoolery, Varian Associates, Palo Alto, Calif., 1962 and 1963.
15. Carbon-13 NMR Spectra. L F Johnson and W C Jankowski, John Wiley and Sons, Inc., New York, N.Y., 1972.
16. Feigl, F: Spottests in organic analysis. 7 Ed. Amsterdam 1966.
17. Bilson, K: A qualitative scheme for the identification of some explosives and related compounds by TLC. MQAD Report No. 212. Feb. 1974.

DATE  
FILMED

9-8